

AERODROME CHART

BEARINGS ARE MAGNETIC.
ALTITUDES, DISTANCES,
ELEVATIONS AND HEIGHTS
IN METERS.

D-ATIS(Chinese) 127.225
D-ATIS(English) 128.4


Delivery01 121.875(121.775)(DCL AVBL)
Delivery02 122.825(DCL AVBL)

TWR01 118.825(124.35)
TWR02 118.375(124.35)
TWR03 130.425(130.3)
TWR04 118.725(130.3)

GND01 121.975(121.775)
GND02 121.625(121.775)
GND03 121.7
GND04 122.6

APN01 122.15(121.775)
APN02 122.7(121.775)

ZBAD BEIJING/Daxing
N39° 30.0'E116° 24.0' ELEV 24.9m



Note: Advanced Surface Movement Guidance and Control Systems(A-SMGCS) in operation.
Mode S transponder and ADS-B(if equipped) is required on at all time when maneuvering on the airport.

RWY	Direction	Bearing strength
17L/R,19R	181°	PCR 1010/R/A/W/T: RWY01L/19R CONC PCR 1000/R/A/W/T: RWY11L/29R CONC PCR 880/R/A/W/T: RWY17L/35R CONC PCR 840/R/A/W/T: RWY17R/35L CONC PCR 1170/R/A/W/T: TWY B9, C(north of M), K10-K14 PCR 1160/R/A/W/T: TWY T3, V PCR 1110/R/A/W/T: TWY B7, T6, T7(east of W6) PCR 1090/R/A/W/T: TWY K3-K5 PCR 1070/R/A/W/T: TWY A10, A12, B6 PCR 1060/R/A/W/T: TWY K(west of K2) PCR 1050/R/A/W/T: TWY C(T5-M) PCR 1040/R/A/W/T: TWY C1, J10, J12, J14, P, P9, Q, Q7-Q9(west of F4), T1, T2, U8, U9 PCR 1030/R/A/W/T: TWY B1, B2, B8, K1, K2, M1-M4 PCR 1010/R/A/W/T: TWY B(south of T5), B3-B5, C(south of T5), E, E0-E9, E12, E13 PCR 1000/R/A/W/T: TWY K(east of K2), L, L2-L7, M(east of L3) PCR 990/R/A/W/T: TWY M0, Z7-Z9 PCR 980/R/A/W/T: TWY C8, J9, J11, J13 PCR 970/R/A/W/T: TWY D, F6, G7-G9, Q9(east of F4), T7(west of W6), T8, T9, W1-W6 PCR 950/R/A/W/T: TWY G, V12-V14, V16, V17 PCR 940/R/A/W/T: TWY M(west of L3) PCR 930/R/A/W/T: TWY B(north of T5), T4, Z1-Z6 PCR 920/R/A/W/T: TWY C2-C4, G3, G5 PCR 900/R/A/W/T: TWY E10, E11, G0-G2, H3-H6, T, Y5(north of Y6), Y6-Y9 PCR 890/R/A/W/T: TWY T5, Y0-Y4, Y5(south of Y6), Z0 PCR 870/R/A/W/T: TWY A, C5-C7 PCR 860/R/A/W/T: TWY G4, G6 PCR 850/R/A/W/T: TWY H, H10, H11 PCR 830/R/A/W/T: TWY J PCR 810/R/A/W/T: TWY A1, A3, A5, A7, J1, J3, J5, J7 PCR 770/R/A/W/T: TWY F4, F5, Q7(east of F4), Q8(F3-F4) PCR 750/R/A/W/T: TWY J2, J4, J6, J8 PCR 740/R/A/W/T: TWY A2, A4, A6, A8 PCR 650/R/A/W/T: TWY F1-F3, Q8(east of F3)
35L/R,01L	001°	
11L	111°	
29R	291°	

Changes: ILS FREQ.