

AERODROME CHART

D-ATIS(English) 128.65
D-ATIS(Chinese) 127.6
TWR 124.3(118.3) TWR01 for 18R/36L
118.5(118.05) TWR02 for 18L/36R
118.6(118.3) TWR03 for 01/19


Delivery01 121.6 (west of RWY 18L/36R) (OCL AVBL)
Delivery02 121.65 (east of RWY 18L/36R) (OCL AVBL)

GND 121.9(121.95) GND01
121.8(121.95) GND02
121.7(121.95) GND03
121.75(121.95) GND04
121.85(121.95) GND05

APN01 122.225 (121.95)
APN02 122.625 (121.95)
APN03 122.675 (121.95)
APN04 122.125 (121.95)

ZBAA BEIJING/Capitol
N40°04.4'E116°35.9' ELEV 35.3m

BEARINGS ARE MAGNETIC.
ALTITUDES, DISTANCES,
ELEVATIONS AND HEIGHTS
IN METERS.



Remarks:

- 1. ARP at center of RWY18L/36R
- 2. Symbol " --- " : HP Holding point
- 3. TWY C4 is used by A/C turn to north FM TWY P4;
- 4. TWY C5 is used by A/C turn to south FM TWY P5;
- 5. --- : GP criticalarea and sensitive area,A/C forbidden to enter without ATC clearance.
- 6. [HS1] : Hot spot.
- 7. x x : U/S TWY.

RWY		Direction	Bearing strength	
01/36L/36R	001°		PCR 930/R/A/W/T: RWY01/19 CONC PCR 1010/R/B/W/T: RWY18L/36R ASPH PCR 1020/R/B/W/T: RWY18R/36L CONC-ASPH PCR 1410/R/B/W/T: TWY D2(S of Z10) PCR 1360/R/B/W/T: TWY Z10 PCR 1340/R/B/W/T: TWY P8, P9 PCR 1300/R/B/W/T: TWY F(BTN Z2 & S4) PCR 1290/R/B/W/T: TWY Z0 PCR 1280/R/B/W/T: TWY C6-C8 PCR 1270/F/B/X/T: TWY P2, P6, P7 PCR 1260/R/B/W/T: TWY P1 PCR 1250/F/B/X/T: TWY A0, A1, A8, A9, E0-E8, W2, W7 PCR 1250/R/B/W/T: TWY Z3(N of Z4) PCR 1240/R/A/W/T: TWY W3-W6 PCR 1240/R/B/W/T: TWY D2(N of Z10), S3 PCR 1230/R/B/W/T: TWY D3-D6, P0 PCR 1220/R/B/W/T: TWY M2-M6 PCR 1200/R/B/W/T: TWY F(N of S4) PCR 1170/R/B/W/T: TWY M PCR 1150/R/B/W/T: TWY M7, Y8, Y9 PCR 1140/R/B/W/T: TWY Z8 PCR 1130/R/B/W/T: TWY Z4, Z6	PCR 1110/R/B/W/T: TWY Y3, Y6 PCR 1060/R/B/W/T: TWY D1, Z7, Z9, Z18 PCR 1050/R/B/W/T: TWY C1, C2, J2, J3, Z3(S of Z2) PCR 1040/R/B/W/T: TWY C PCR 1030/R/B/W/T: TWY D7, D8, J(N of U2), S4, S5 PCR 1020/R/A/W/T: TWY M0, M1 PCR 1020/R/B/W/T: TWY Q2-Q7, Z3(BTN Z2 & Z4) PCR 1010/R/B/W/T: TWY F(S of Z2), Z2 PCR 960/R/B/W/T: TWY Z12 PCR 940/R/A/W/T: TWY Q0, Q1, Q8, Q9, U2 PCR 920/R/B/W/T: TWY C3-C5 PCR 910/R/A/W/T: TWY J(S of U2), J1, J4, K, K0-K7, S6, S7 PCR 900/R/A/W/T: TWY U3-U9 PCR 900/R/B/W/T: TWY P3 PCR 880/R/A/W/T: TWY F0-F4, F7, J5, J6, T1-T6 PCR 870/R/A/W/T: TWY W0, W9, Y1, Y2, Y4, Y5, Y7 PCR 850/R/A/W/T: TWY G, G0-G7, H, H0-H2, H4-H7 PCR 820/F/B/X/T: TWY P4, P5 PCR 780/R/B/W/T: TWY Z15, Z16 PCR 710/R/A/W/T: TWY Z1 PCR 590/R/B/W/T: TWY Z20-Z24

Changes: PCR, APRON NI.