

AERODROME CHART

D-ATIS(English) 127.85
D-ATIS(Chinese) 128.65


TWR01 118.8(118.325) for RWY17L/35R, 17R/35L
TWR02 118.4(118.725) for RWY16L/34R, 16R/34L
TWR03 124.35(118.325) for RWY17L/35R
TWR04 118.575(118.725) for RWY16R/34L

APN01 121.65(122.125)
APN02 121.975(122.125)
APN03 122.7(122.125)
APN04 122.6(122.125)

GND01 121.7
GND02 121.8
GND03 121.875
GND04 121.625
Delivery: 121.95(121.625)(DCL

ZSPD SHANGHAI/Pudong

N31°08.7'E121°47.6' ELEV 3.8m



Note:
RWY 16R/34L MARKING AND LIGHTING AIDS
SIMILAR TO RWY 17L/35R.
RWY 16L/34R MARKING AND LIGHTING AIDS
SIMILAR TO RWY 17R/35L.
ARP at center of RWY 17L/35R.

RWY	Direction	RWY and TWY bearing strength		
16L/ 16R/ 17L/ 17R	168°	PCR 980/R/B/W/T: RWY16L/34R CONC PCR 960/R/B/W/T: RWY16R/34L CONC PCR 890/R/A/W/T: RWY17L/35R CONC PCR 930/R/A/W/T: RWY17R/35L CONC	PCR 1780/F/C/X/T: TWY P1(BTN TWY D & RWY17L/35R), P2-P6(BTN TWY A&D) PCR 1680/F/C/X/T: TWY R1(BTN TWY G & RWY16R/34L), R2(BTN TWY G & RWY16R/34L), R3(E of TWY F), R4(E of TWY F), R5(BTN TWY G & RWY16R/34L), R6(BTN TWY G & RWY16R/34L) PCR 1620/F/C/X/T: TWY G(N of TWY G6 exit & S of TWY G1 exit) PCR 1100/R/A/W/T: TWY B7 PCR 1030/R/B/W/T: TWY S1, S2, T5, T6 PCR 1000/R/A/W/T: TWY A1-A4 PCR 980/R/A/W/T: TWY B2-B6, F1-F4 PCR 960/R/A/W/T: TWY B8 PCR 950/R/A/W/T: TWY B1, C1, C2, C5, C6, D1, D2, D5, D6 PCR 940/R/A/W/T: TWY A5, A6 PCR 910/R/A/W/T: TWY A, B PCR 900/R/A/W/T: TWY P1(BTN TWY B & RWY17L/35R) PCR 890/R/A/W/T: TWY H1-H6, J1, J2 PCR 880/R/A/W/T: TWY Q1-Q6	PCR 870/R/A/W/T: TWY E0-E8, W1-W7 PCR 860/R/A/W/T: TWY G2-G5 PCR 850/R/A/W/T: TWY G1, G6 PCR 840/R/A/W/T: TWY P3(E of TWY A), P4(E of TWY A & W of TWY D) PCR 820/R/A/W/T: TWY H, P2(E of TWY A & W of TWY D), P5(E of TWY A) PCR 810/R/A/W/T: TWY C3, C4, D3, D4, R2(E of TWY G & W of RWY16R/34L), R3(W of TWY F), R4(W of TWY F), R5(E of TWY G & W of RWY16R/34L) PCR 790/R/A/W/T: TWY C, D, E, F, T3, T4, V3-V8 PCR 780/R/A/W/T: TWY G(BTN TWY G1&G6 exit) PCR 770/R/A/W/T: TWY P6(E of TWY A & W of TWY D) PCR 750/R/A/W/T: TWY R1(E of TWY G & W of RWY16R/34L), R6(E of TWY G & W of RWY16R/34L) PCR 740/R/A/W/T: TWY P1(W of TWY D), T1, T2, V1, V2
34R/ 34L/ 35R/ 35L	348°			

Changes: Nil.